MicroPort Minos™ Sterile

Abdominal Aortic Stent-Graft and Delivery System

INSTRUCTIONS FOR USE

EN	<i>Abdominal Aortic Stent-Graft and Delivery System</i>	1
ES	<i>Sistema de injerto y colocación de endoprótesis aórtico abdominal</i>	16
IT	<i>Sistema di endoprotesi aortico addominale e sistema di rilascio</i>	33
DE	<i>Abdominales Aortenstentgraft- und Einführsystem</i>	50
CS	<i>Stentgraft pro břišní aortu a zaváděcí systém</i>	68
SK	<i>Abdominálny aortálny stentgraft a zavádzací systém</i>	83
FR	<i>Système de pose du Stent-Graft pour l'aorte abdominale</i>	98
HU	<i>Hasi aorta sztentgraft és felhelyező rendszer</i>	115
BG	<i>Стент-графт за абдоминална аорта и система за въвеждане</i>	130
RO	<i>Sistem de greșă de tip stent aortic abdominal și de plasare</i>	148

[SYMBOL DESCRIPTION]

1)		Do not re-use;No re utilizar;Non riutilizzare;Nicht wiederverwenden;Nepoužívejte opakovaně;Nepoužívajte opakovane;Ne pas réutiliser;Újrahasználni tilos!;Да не се използва повторно;A nu se reutiliza
2)		Do not re-sterilize;No re esterilizar;Non risterilizzare;Nicht erneut sterilisieren;Nesterilizujte opakovaně;Nesterilizujte opakovane;Ne pas restériliser;Újrasterilizálni tilos!;Не стерилизирайте повторно;Nu reesterilizați
3)		Consult instructions for use;Consulte las instrucciones de uso;Consultare le istruzioni per l'uso;Anwendungshinweise beachten;Prostudujte si pokyny k použití;Prečítajte si návod na použitie;Lire la notice d'utilisation;Olvasa el a használati utasítást;Вижте инструкциите за употреба;Consultați instrucțiunile de utilizare
4)		Keep away from sunlight;Mantener alejado de la luz solar;Tenere lontano dalla luce solare;Vor Sonneneinstrahlung schützen;Udržujte mimo dosah slunečního světla;Chráněte pred slnečným žiarením;Conserver à l'abris de la chaleur;Napfénytől távol tartandó;Пазете от слънчева светлина;A se feri de razele soarelui
5)		Keep dry;Mantener seco;Mantenere asciutto;Trocken halten;Udržujte v suchu;Uchovávejte v suchu;Maintenir au sec;Szárason tartandó;Пазете сух;A se păstra în stare uscată
6)		Use-by date;Utilizar por fecha;Da utilizzare entro il;Verfallsdatum;Použit do data;Dátum spotreby;Date limite d'utilisation;Felhasználható eddig a dátumig;Срок на годност;Valabil până la
7)		Date of manufacturer;Fecha de fabricación;Data di fabbricazione;Herstellungsdatum;Datum výroby;Dátum výrobcu;Date de fabrication;Gyártás időpontja;Дата на производство;Data fabricantului
8)		Batch code;Código de lote;Codice lotto;Chargennummer;Kód šarže;Kód šarže;Numéro de lot;Gyártási tételszám;Код на партидата;Cod de lot
9)		Catalog number;Numero de catalogo;Numero catalogo;Katalognummer;Katalogové číslo;Katalógové číslo;Numéro de catalogue;Katalógusszám;Каталожен номер;Număr de catalog
10)		Sterilized using ethylene oxide;Esterilizado con óxido de etileno;Sterilizzato mediante ossido di etilene;Sterilisiert mit Ethylenoxid;Sterilizováno pomocí etylenoxidu;Sterilizované etylénoxidom;Stérilisé à l'oxyde d'éthylène;Etilén-oxiddal sterilizálva;Стерилизирано с етиленов оксид;Sterilizat cu oxid de etilenă
11)		Do not use if package is damaged;No lo use si el paquete está dañado;Non utilizzare se la confezione è danneggiata;Bei beschädigter Verpackung nicht verwenden;Nepoužívejte, pokud je obal poškozen;Nepoužívajte, ak je obal poškodený;Ne pas utiliser si le conditionnement est endommagé;Ne használja, ha sérült a csomagolás;Не използвайте, ако опаковката е повредена.;A nu se utiliza în cazul în care ambalajul este deteriorat
12)		Contents: One Unit;Contenido: Una unidad;Contenuto: Una unità;Inhalt: Eine Einheit;Obsah: 1;Obsah: 1;Contenu : 1;Tartalma: 1;Съдържание: Един елемент;Conținut: O unitate
13)		Manufacturer;Fabricante;Produttore;Hersteller;Výrobce;Výrobca;Fabricant;Gyártó;Производител;Producător

14)	<div style="border: 1px solid black; padding: 2px; display: inline-block;"> EC REP </div>	<p>Authorized representative in the European Community; Representante autorizado en la Comunidad Europea; Rappresentante autorizzato nella Comunità Europea; Autorisierter Vertreter für die Europäische Gemeinschaft; Zplnomocněný zástupce v Evropském společenství; Splnomocnený zástupca v Európskom spoločenstve; Représentant autorisé dans la Communauté Européenne; Hivatalos képviselő az Európai Közösségben; Упълномощен представител в Европейската общност; Reprezentant autorizat în Comunitatea Europeană</p>
-----	--	--

[PRODUCT DESCRIPTION]

MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System (hereinafter referred to as the graft system) is for the purpose of interventional treatment of abdominal aortic aneurysm.

Abdominal aortic aneurysm is shown in Figure 1.


Figure 1 Schematic Diagram of Abdominal Aortic Aneurysm

MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System treats abdominal aortic aneurysms using microtrauma techniques and principle of intracavitary isolation.

Endovascular Stent Graft is separately structured, consisting of a main body and limb stent grafts. They are successively delivered into the lesion from the bilateral femoral artery through corresponding delivery system, and then released and connected as a whole in-vivo to form a new blood flow channel and exclude the aneurysm from the blood circulation system; thereby the aneurysm cavity is isolated, and the purpose of treating the aneurysm is achieved.

If necessary, a CUFF stent graft can be placed at the proximal end of the main body stent graft or an additional limb stent graft can be placed at the distal end of the limb stent graft as an extension.


The Endovascular Stent Graft is self-expanding graft and made of nitinol, which has super-elasticity and can support the stent graft strongly. The stent's material is a linear polyester textile (commonly known as polyester), which is a general-purpose artificial blood vessel material with good biocompatibility.

[SPECIFICATION AND MODEL]


When determining the size specifications of the Endovascular Stent Graft, the relevant dimensions of the abdominal aortic aneurysm should be accurately measured, such as: the length and diameter of the proximal end of aneurysm neck, the length of the infrarenal abdominal aorta, the length of the bilateral common iliac artery, the range of effect of the aneurysm of the common iliac artery, and the diameter of the distal end of aneurysm neck, etc.

MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System's models are indicated by letters and are classified into three types: the model of the main body stent graft system as CM; the model of the CUFF stent graft system as CC; and the model of the limb stent graft system as CL. The specifications corresponding to each model are indicated by diameters D1, D2, D3 and lengths L1, L2, and L3 of the


respective stent grafts. D1 represents the proximal end diameter of the main body stent graft, D2 represents the diameter of the CUFF stent graft, D3 represents the diameter of the distal end of the stent graft; L1 represents the ipsilateral stent length of the main body stent graft, L2 represents the stent length of the CUFF stent graft, and L3 represents the stent length of the limb stent graft.

Main body stent graft

CUFF stent graft

Limb stent graft

Figure 2 Schematic Diagram of the Stent Graft

Combined with the schematic diagram of the stent graft in Figure 2, the product model specification code is exemplified as follows.

1) CM22-90

“CM” stands for the main body stent graft system

“22” stands for the proximal end diameter of the main body stent graft (D1)

“90” stands for the ipsilateral stent length of the main body stent graft (L1)

***The contralateral length of the main body stent graft is fixed at 80 mm, and the diameter of the ipsilateral and contralateral distal ends is fixed at 12mm.**

2) CC22-40

“CC” stands for CUFF stent graft system

“22” stands for the diameter of the CUFF stent graft (D2)

“40” stands for the length of the CUFF stent graft (L2)

3) CL13-80

“CL” stands for the limb stent graft system

“13” stands for the distal end diameter of the limb stent graft (D3)

“80” stands for the stent length of the limb stent graft (L3)

***The diameter of the proximal end of the stent graft is fixed at 14 mm**

The model specifications for MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System and the mating relationship between the graft and the delivery system specifications are shown in Table 1.

Table 1 Model and specifications of MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System

Model and specifications of main body stent graft system	Profile of the delivery system	Working length of the delivery system
CM22-90	14F	550 mm
CM22-100		
CM22-110		
CM22-120		
CM22-130		
CM22-140		
CM24-90		
CM24-100		
CM24-110		
CM24-120		
CM24-130		
CM24-140		
CM26-90		
CM26-100		
CM26-110		
CM26-120		
CM26-130		
CM26-140		

Model and specifications of main body stent graft system	Profile of the delivery system	Working length of the delivery system
CM28-90	16F	550 mm
CM28-100		
CM28-110		
CM28-120		
CM28-130		
CM28-140		
CM30-90		
CM30-100		
CM30-110		
CM30-120		
CM30-130		
CM30-140		
CM32-90		
CM32-100		
CM32-110		
CM32-120		
CM32-130		
CM32-140		
CM34-90		
CM34-100		
CM34-110		
CM34-120		
CM34-130		
CM34-140		

Model and specifications of CUFF stent graft system	Outer diameter of the delivery system	Effective length of the delivery system
CC22-40	16F	550mm
CC24-40		
CC26-40		
CC28-40		
CC30-40		
CC32-40		
CC34-40		
CC36-40		

Model and specifications of limb stent graft system	Outer diameter of the delivery system	Effective length of the delivery system
CL10-80	12F	610mm
CL10-90		
CL10-100		
CL10-110		
CL10-120		
CL10-130		
CL10-140		
CL13-80		
CL13-90		
CL13-100		

Model and specifications of limb stent graft system	Outer diameter of the delivery system	Effective length of the delivery system
CL13-110	12F	610mm
CL13-120		
CL13-130		
CL13-140		
CL16-80		
CL16-90		
CL16-100		
CL16-110		
CL16-120		
CL16-130		
CL16-140		
CL18-80		
CL18-90		
CL18-100		
CL18-110		
CL18-120		
CL18-130		
CL18-140		
CL20-80		
CL20-90		
CL20-100		
CL20-110		
CL20-120		
CL20-130		
CL20-140		
CL22-80		
CL22-90		
CL22-100		
CL22-110		
CL22-120		
CL22-130		
CL22-140		
CL24-80		
CL24-90		
CL24-100		
CL24-110		
CL24-120		
CL24-130		
CL24-140		

The delivery systems of the main body stent graft, limb stent graft and CUFF stent graft are basically structured similarly. The effective length of the outer sheath of the delivery system of the main body stent grafts and CUFF stent grafts is 55cm, the effective length of the outer sheath of the limb stent grafts delivery system is 61cm, and the specifications of outer sheath of the main body stent grafts delivery system consist of 14F and 16F. The specification of outer sheath of the CUFF stent grafts delivery system is 16F, and the specification of outer sheath of the limb stent grafts delivery system is 12F. Prior to leaving the factory, the stent grafts have been individually assembled in the delivery systems and

separately packaged and sterilized. The delivery system is to deliver the stent graft to the lesion and release it to appropriate site by operating the delivery system. The structure of the graft system is shown in Figure 3.

[INDICATIONS]

MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System is applicable to abdominal aortic aneurysm with proximal aneurysm neck length ≥ 15 mm.

[CONTRAINDICATIONS]

- There is severe stenosis or calcification near the aneurysm, which may easily lead to the difficulty in adhering the stent graft;
- There is severe stenosis, distortion at the iliac artery or femoral artery that cannot provide effective access to the product;
- Patients with ruptured aneurysms or acute aneurysms;
- Patients with infectious aneurysms;
- Those who are allergic to contrast agents or renal insufficient that cannot tolerate contrast agents;
- Patients with severe coagulopathy and is prone to increase post-operative bleeding;
- The aneurysm neck bending angle is $>60^\circ$; no enough anchoring area at the distal end, or those who are prone to cause the stent graft hardly adherent;
- The blood supply branch of the vital organ of the patient starts from the aneurysm cavity;
- The patient has other diseases and the life expectancy is less than one year;
- Suffering from connective tissue disease, such as a patient with Marfan syndrome;
- Minors and pregnant women;
- Those who are allergic to nitinol;
- Others who are inapplicable for stenting.

[WARNING]

- Please read the Instructions for Use (IFU) carefully before use to avoid adverse consequences.
- Before using, physicians need to undergo appropriate professional training to have a clear understanding about the principles, clinical applications, complications, side effects and damage of aortic stenting;
- The product consists of a main body stent graft and delivery system, a limb stent graft and delivery system, and a CUFF stent graft and delivery system, each packaged separately;
- MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System has been sterilized by ethylene oxide before leaving the factory;
- The product is disposable. The company takes no responsibility for any cross-contagion or other adverse consequences resulted from repeated use of the stent graft, delivery system or any other components therein.
- Do not use if the package is damaged before unpacking.
- The user shall not resterilize the product after opening the package.


Figure 3 Schematic diagram of the structure of the abdominal major stent graft system

[COMPLICATION]

The use of a stent graft and delivery system may result in the following or other complications:

- Endoleak;
- Aneurysm rupture;
- Iliac artery or femoral artery injury;
- Distal vascular embolization of the limb;
- Graft thrombosis;
- Intestinal ischemic necrosis;
- Sacrococcygeal ischemic necrosis;
- Renal obstruction;
- Wound infection;
- Graft infection;
- Arteriovenous fistula;
- Lymphorrhagia;
- Paraplegia.

[SPECIFICATIONS SELECTION]

It is critical to select the specifications of the stent graft that improper selection of the graft diameter or stent segment length can result in surgical failure. When determining the specifications, the diameter and length of the corresponding parts in Figure 1 should be accurately measured. The diameter D1 of the main body portion of the main body stent graft should be about 15% larger than the diameter of the proximal end/distal end neck of the aneurysm measured in Figure 1. It is recommended to refer to Table 2 for selection. The total length of the main body stent graft L1 should ensure effective coverage over the abdominal aortic aneurysm, but not more than the opening distal end of the lower of the two renal arteries; the total length of the limb stent graft L3 should ensure effective coverage over the aneurysm after being nested with the main body stent graft (see Figure 2).

Table 2 Suggested correspondence relationship between graft diameter and target vessel inner diameter (ID)

Main body stent graft proximal end diameter D1(mm)	Suggested target vessel ID (mm)
22	17-19
24	20-21
26	22-23
28	24-25
30	26-27
32	28-29
34	30-31

CUFF stent graft diameter D2 (mm)	Suggested target vessel ID (mm)
22	17-19
24	20-21
26	22-23
28	24-25
30	26-27
32	28-29

CUFF stent graft diameter D2 (mm)	Suggested target vessel ID (mm)
34	30-31
36	32-33

Diameter of the distal end of the limb stent graft D3(mm)	Suggested target vessel ID (mm)
10	8-9
13	10-12
16	13-14
18	15-16
20	17-18
22	19-20
24	21-22

[ACCESSORIES REQUIRED]

- 1) Free angle C-arm and equipped with:
High resolution fluoroscopy equipment;
High quality angiography equipment;
Digital subtraction angiography equipment (DSA);
- 2) Surgical equipment to respond to unexpected situations;
- 3) Tweezers;
- 4) Wire. It is recommended to use a 0.035" J-shaped super-hard guidewire with a length ≥ 2.6 m;
- 5) Percutaneous needles and vasodilator devices;
- 6) Balloon catheter;
- 7) Angiographic catheter;
- 8) Contrast agent;
- 9) Heparin saline.

[PREOPERATIVE PREPARATION]

Adjust the C-arm of the bed device to proper position; in vitro, observe the graft installation area of the main body graft system through the monitor, rotate the delivery system to confirm the orientation of the lateral radiopaque marker on the contralateral branch of the main body stent graft and keep it contralateral the patient during subsequent procedures. When the main body stent graft is installing, the ipsilateral branch is always on the same side as the infusion tube of the delivery system (refer to Figure 3).

[SURGICAL STEPS]

- 1) Anesthetize the patient.
- 2) Puncture or incise in the patient's inguinal region.

Note: *in principle, the access route of the main body graft delivery system should be determined in compliance with the following three points: 1, the diameter of the lateral common iliac artery / external iliac artery, femoral artery is relatively large; 2, the lateral common iliac artery/ external artery and femoral artery are smooth and no severe bending; 3, the common iliac artery and the abdominal aorta angulate small.*

- 3) Deliver a J-type super-hard guidewire with a diameter of 0.035"

until it is above the opening of the celiac trunk artery.

- 4) Contrast. Place the contrast catheter at a position above the double renal artery for angiography, analyze the lesion, mark the renal artery opening and other important parts, and confirm the size and specifications of the stent graft.

Note:after angiography, keep the patient's position not moving until the main body stent graft is released.

Note:if necessary, use a graduated contrast catheter to review the size and specifications of the stent graft.

- 5) Graft system preparation.
 - a. Remove the main body graft system of the MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System of the predetermined specification from the package;
 - b. Inspect the product for any damage; replace the product if there is damage and continue to use if the product is normal;
 - c. Tilt the delivery system so that the tapered head is obliquely pointed upward; inject heparin saline into the delivery system from the Exchange Sheath Flushing Port to fill the entire outer sheath;
 - d. Close the Exchange Sheath Flushing Port of the delivery system;
 - e. Process the branch graft system of the predetermined specification in the same manner as described above.
- 6) Deliver in the main body graft system. Insert the super-hard guidewire into the tapered head of the main body stent graft system, hold the outer sheath and push the system along the super-hard guidewire into the appropriate position in the patient.

Note:the infusion tube should be at the lateral of the patient during introducing into the blood vessel.

- 7) Determine the position of the main body stent graft. Use the contrast agent, the marking line on the monitor, or the scale on the patient's surface to ensure the following positions:
 - a. When the lower edge of the proximal end radiopaque marker of the stent graft is at the distal end of the lowermost renal artery opening, the most distal radiopaque marker on the contralateral branch is higher than the abdominal aortic bifurcation, but not in the ipsilateral common iliac artery;
 - b. The radiopaque marker of the contralateral branch is in the direction of the contralateral radial artery.

Note:this process is critical, otherwise it may lead to surgery failure!

- 8) Release the main body stent graft. Hold the proximal handle, push the handle switch and rotate the rear handle in the direction of the handle arrow, and release the contralateral branch of the main body stent graft. After the contralateral branch of the main body graft is ejected from the outer sheath, stop rotating, so that the ipsilateral branch of the main body graft remain in the outer sheath; release the switch after rotating, pull the rear release puller to release the proximal end of the graft.

Note:this step must be carried out with caution in the

dynamic image, otherwise it may lead to surgery failure!

Note:excessive blood pressure can increase the difficulty for positioning the stent graft, which may result in inaccurate positioning of the graft during release, thus before starting to release, ensure that the patient's blood pressure is controlled at around 80 mmHg.

Note:minimize the bending of the exposed human part of the delivery system, otherwise it may result in inaccurate positioning of the graft during release.

Note:during the release, avoid unintentional movement or rotation of the delivery system. Do not pull the rear release puller before the contralateral branch is completely released from the outer sheath.

- 9) Establish a contralateral branch guidewire channel. Introduce the super-hard guidewire from the contralateral into the main body stent graft through the contralateral branch of the main body stent graft. The optional ways to enter the guidewire are as follows:
 - a. Directly send the guidewire from the contralateral side into the contralateral branch until it exceeds the renal artery opening;
 - b. Puncture from the left subclavian artery and deliver the exchange guidewire and catheter into the contralateral branch and lead out from the contralateral incision. After withdrawing the exchange guidewire, deliver the super-hard wire along the catheter.

Note:make sure that the super-hard guidewire passes through the inside of the main body stent graft, not between the stent graft and the aneurysm wall, otherwise it will lead to surgical failure. It can be confirmed by delivering the contrast catheter from the contralateral side and performing angiography.

Note:throughout this step, make sure that the main body delivery system is always fixed.

- 10) Release the ipsilateral branch of the main body stent graft. Hold the proximal handle of the delivery system, rotate the rear handle in the direction of handle arrow, and release the remaining part of the ipsilateral side of the main body stent graft. Rotate the rear handle to fully release the main body graft.

Note:this step must be carried out with caution in the dynamic image, otherwise it may lead to surgery failure!

Note:the outer sheath marker band should be 20 mm or more away from the distal end of the graft to ensure that the graft has been fully released.

- 11) Withdraw the main body graft system. Keep the super-hard guidewire still, press the quick connector switch, hold the proximal end of the quick connector (front outer sheath) still, move

the distal end of the quick connector (rear outer sheath). Withdraw the tapered head of the delivery system from the body along the super-hard wire, leaving the front outer sheath in the body.

Note:*The super-hard guidewire should remain in the body after the delivery system has been withdrawn.*

- 12) Introduce the contralateral branch graft system. Introduce the branch graft system along the contralateral super-hard wire into the contralateral branch of the main body stent graft; adjust the position of the branch graft system to make sure the distal radiopaque marker of the branch graft is in the desired position.

Note:*after positioning, ensure that the limb stent graft is inside the contralateral branch of the main body stent graft, and the proximal end radiopaque marker of the limb stent graft is not higher than the radiopaque marker of the bifurcation of the stent graft, not lower than the second radiopaque marker at the proximal end of the contralateral branch of the main body.*

- 13) Release the contralateral branch graft. Hold the proximal handle, push the handle switch and rotate the rear handle in the direction of handle arrow to release the graft completely from the outer sheath. After the graft is ejected from the outer sheath, rotate the rear release switch. Pull the rear release puller to extreme position and release the proximal end of the branch graft.

Note:*this step must be carried out with caution in the dynamic image, otherwise it may lead to surgery failure!*

Note:*minimize the bending of the exposed human part of the delivery system, otherwise it may result in inaccurate positioning of the graft during release.*

Note:*during the release, avoid unintentional movement or rotation of the delivery system. Do not pull the rear release puller before the contralateral branch is completely released from the outer sheath.*

- 14) Withdraw the branch graft system. Keep the super-hard wire still, press the quick connector switch, hold the proximal end of the quick connector (front outer sheath) still, move the distal end of the quick connector (rear outer sheath). Withdraw the tapered head of the delivery system from the body along the super-hard guidewire, leaving the front outer sheath in the body

Note:*the super-hard guidewire should remain in the body after the delivery system has been withdrawn.*

- 15) Deliver along the ipsilateral super-hard guidewire into the ipsilateral branch graft system; dock the front outer sheath of the branch system on the quick connector of the main body system and then enter the ipsilateral branch of the main body stent graft; adjust the position of the graft system to make sure the distal radiopaque marker of the branch graft is in the desired position.

Note:*make sure the limb stent graft is inside the main*

stent graft branch.

Note:the proximal end radiopaque marker of the limb stent graft is not higher than the radiopaque marker of the bifurcation of the stent graft, and not lower than the first radiopaque marker at the proximal end of the ipsilateral branch of the main body.

Note:during the introduction of the ipsilateral branch graft system, avoid moving the front outer sheath of the main body system which may cause displacement to the distal end of the main body graft.

- 16) Release the ipsilateral limb stent graft. Hold the proximal handle, push the handle switch and rotate the rear handle in the direction of handle arrow to release the graft completely from the outer sheath. After the graft is ejected from the outer sheath, rotate the rear release switch. Pull the rear release puller to extreme position and release the proximal end of the branch graft.

Note:this step must be carried out with caution in the dynamic image, otherwise it may lead to surgery failure!

- 17) Withdraw the ipsilateral limb stent graft delivery system. Keep the super-hard guidewire still, after the tapered head is retracted into the outer sheath of the branch system along the super-hard guidewire, withdraw the delivery system completely from the body.

- 18) Evaluate the effect of graft release. Perform angiography to evaluate the shape, position and closure effect of the stent graft after release. If the effect is not satisfactory, the below follow-up measures should be taken:

- a. If the proximal end position of the stent graft is suitable, but there is obvious type I endoleak, deliver a balloon of corresponding size in time to properly expand the proximal end stent segment of the graft and to make the proximal end stent fully adherent;
- b. If the proximal end position of the stent graft is lower, and the endoleak is not eliminated after using a balloon, or the proximal end stent graft fails to effectively isolate the aneurysm cavity, a suitable CUFF stent graft, whose release operation is consistent with the release of the main body graft, should be released at the proximal end of the main body stent graft.
- c. If the position of the limb stent graft is appropriate, but there is an obvious type III endoleak at the overlap with the main body, a balloon of corresponding size should be fed to expand the overlapping area appropriately to make the stents of the main body and the branch fully fit;
- d. If the distal position of the limb stent graft is too high and not effectively isolate the aneurysm cavity, add a suitable extension to the distal end of the limb stent graft.

Note:please operate with caution when taking follow-up actions to avoid displacement of the stent graft.

Note:after completing the angiography, a guidewire should be kept in the body for subsequent operation.

19) Withdraw all the guidewires, suture the incision, and complete the surgery.

[MRI INFORMATION]

Non-clinical testing indicates that MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System is safe in the following MRI (Magnetic Resonance Imaging) environments:

Philips Achieva 3.0T TX, magnetic field strength 3-Tesla, cavity diameter=60cm, active shielding, short magnet (magnet length 1.57m), maximum gradient field strength 80mT/m, maximum gradient switching rate 200mT/m/ms, software: Philips DICOM Viewer, version: R2.6L1-SP1;

Magnetically induced displacement force: the maximum deflection angle is 6.5°, less than 45°, that is, the maximum magnetic displacement force is less than gravity;

Radio frequency induced heating: WSAR=3.7W/kg, the temperature rise caused by stent graft upon 20 minutes scanning is 1.73°C;

Magnetically induced torque: no torque;

Image artifacts: Image quality may be affected by MRI-evaluated development area when it is located at and near the graft.

[HOW SUPPLIED]

STERILE: FOR SINGLE USE ONLY.

CONTENTS

- One (1) MicroPort Minos™ Abdominal Aortic Stent-Graft and Delivery System;
- One (1) Instructions for Use.

STORAGE AND TRANSPORTATION

The product should be stored in a dry, clean, well ventilated place under room temperature and noncorrosive atmosphere. Protect the product from heavy pressure, direct sunlight and rain and snow during transportation.

DISPOSAL INSTRUCTIONS

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

[SHELF LIFE]

The shelf life is 2 years under the specified storage conditions. The dates of manufacturer and use-by date can be found on the packaging label.

[DISCLAIMER]

The potential adverse events of the Product listed in the Product Specification are learnt from the information and data collected by Shanghai MicroPort Endovascular MedTech (Group) Co., Ltd. (hereinafter referred to as “Endovastec™”) from the existing clinical trials within People’s Republic of China and other countries and regions, and from the information and data collected after the Product was launched in the market and shall not be deemed as factually occurred. The aforesaid

potential adverse events are features of the Product for which sufficient warning messages have been provided and shall be deemed neither as product defects nor as an exclusive list of all potential risks of the Product. The Product may have other potential risks, which have not been discovered or known due to the current scientific and technological level and the limitation of Endovastec™'s cognitive capability. Such potential risks shall not be deemed as product defects as well.

Endovastec™ shall not be responsible for any medical cost, any direct, incidental or indirect loss or any claim for damages based on warranty, contract, tort or other legal principles that are caused by human error accidents including without limitation misuse (including reuse), improper choice of product specification and erroneous operation, illegally purchased Products, risks that are undiscoverable by current technological capability and adverse events or risks listed in this Product Specification.

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MicroPort Medical B.V.